

# Just Breathe: A Breathing Control App to Manage Well-Being

## *Authors:*

Francesco Bruno  
Gaetano Sferrazza  
Lorenzo Tonelli  
Nicolò Picchi



# OUTLINE

---


- INTRODUCTION
- JUST BREATHE – GUI DESCRIPTION
- ARCHITECTURE
- BREATHING SIGNAL SCORE
- EXPERIMENTAL RESULTS
- CONCLUSION
- VIDEO DEMONSTRATION

# INTRODUCTION

---




- Breathing is an automatic function of the body that is controlled by the respiratory centre of the brain
- Breathing control is a fundamental operation to manage emotions or states of confusion.
- Just Breathe is an intuitive application to follow a breathing pattern
- Just Breathe Objective: help a user to improve its breathing control




# JUST BREATHE – GUI DESCRIPTION

- We recommend Calibration
- Click on info button to have more info about that pattern
- Chose one of three pattern by clicking on one picture
- We suggest to use a watch to take the time of breathing.



# ARCHITECTURE



**1. Recommended calibration**

**2. Pattern choice by graphical interface**

**3. Accelerometer signals processing:**

1. Calculation of the window average
2. Signal filtering through the Butterworth filter



**4. Computation of the respiratory signal score in according to the pattern**

1. Finding of peak function
2. Application of the formula for scoring



**5. Return of the score to the GUI**



# PATTERN 4-0-4 – SIGNAL PROCESSING



# PATTERN 4-7-8 - SIGNAL PROCESSING



## BREATHING SIGNAL SCORE – DESCRIPTION

---

To measure the **goodness** of a user's breathing, a heuristic algorithm has been developed, **breathingScore()**.

- It takes as input the processed breath signal and the pattern to follow.
- Returns a value between 0 and 1 which will be multiplied by 100 before showing it on the smartphone screen.

$$S = \alpha_1 q_{max} + \alpha_2 q_{min} + \alpha_3 q_{int} + \alpha_4 q_{hold}$$

- ***q<sub>max</sub>***: is the index that measures how close the local maximum peaks are to their average
- ***q<sub>min</sub>***: is the index that measures how close the local minimum peaks are to their average
- ***q<sub>int</sub>***: is the index that measures how much the period of inhalation and exhalation corresponds
- ***q<sub>hold</sub>***: is the index that measures the general goodness of the apnea periods of the entire session




## BREATHING SIGNAL SCORE – EXAMPLE

An example of `breathingScore()` that takes as input the signal in the figure with following parameters:

- `peak_max_tolerance` = 0.005 m/s<sup>2</sup>
- `peak_min_tolerance` = 0.005 m/s<sup>2</sup>
- `breath_period_tolerance` = 500ms
- $\alpha_1 = \alpha_2 = 0.25, \alpha_3 = 0.50$

Computing:

- $q_{\min} = 8/8 = 1$
- $q_{\max} = 7/8 = 0.875$
- $q_{\text{int}} = 15/15 = 1$
- $S = 0.25 q_{\min} + 0.25 q_{\max} + 0.50 q_{\text{int}}$   
 $= 0.954 \rightarrow \mathbf{95.4}$



## EXPERIMENT RESULTS – CHOICE OF HYPERPARAMETERS

---

Peak Tolerance	Breath Period Tolerance	$\alpha_1 = \alpha_2$	$\alpha_3$
0.15	1000	0.3	0.4
0.005	500	0.25	0.5

### Useful Considerations:

- Each user reported difficulties in respecting the chosen breathing pattern for over 40 seconds.
- Each user has a different abdominal extensibility and sensitivity
- Two users stated that keeping the breathing rhythm accurate to the ideal one was difficult.

## EXPERIMENT RESULTS – TEST-SET

---


The test-set created to conduct the experiments was defined as follows:

- **Breathing Pattern:** 4-0-4 for all registration sessions
- **Environment Set:** User lying on bed with smartphone resting on abdomen
- **Breathing Session Duration:** At least 45s
- **Session Type:** 50% Correct - 50% Incorrect (compared to the chosen pattern)
- **Beta-Version** of the JustBreathe App to experimentation
- 60 breathing sessions were collected on a sample of 4 male persons, same build and breathing capacity who installed the **JustBreathe** app on their Android Smartphone

## EXPERIMENT RESULTS – MEASUREMENT OF SCORE RELIABILITY

### Sessions Labels:

- “Good”: When the user breathed in the manner closest to the chosen breathing pattern
- “Bad”: When the user did not breathe in the manner close to the chosen breathing pattern



	Precision	Recall	F1-Score	Support
“Good”	0.90	0.87	0.88	30
“Bad”	0.87	0.90	0.89	30
Accuracy			0.88	60

Table: Metrics Results

# CONCLUSION

---



- **Advantages:**

- Original and innovative solution in the mobile app industry dedicated to supporting user breathing
- Reduced complexity due to the fact that no AI technique is required
- Easy to use and based on an unobtrusive approach for the monitoring phase

- **Limitations:**

- JustBreathe App is not yet ready for large-scale use
- Choice of Hyperparameters based on a too small sample of users
- Physiological characteristics of the users
- Noise in the raw signal obtained from the accelerometer

# VIDEO DEMONSTRATION

---







THANKS FOR YOUR  
ATTENTION!

## REFERENCES

---



- [1]. Francesco Bruno; Gaetano Sferrazza; Lorenzo Tonelli; Nicolò Picchi,  
Manage Well-Being Through a Breathing Control App  
[https://github.com/francescoB1997/Just\\_Breathe](https://github.com/francescoB1997/Just_Breathe)